

jasco



BIP 1, 880/890 SERIES

## optimax check valves

optimize P/N	Optimize Description	JASCO P/N
For BIP 1		
42-56-01038	OPTI-MAX Inlet Housing & Cartridge	6260-H101A
42-56-01043	OPTI-MAX Outlet Housing & Cartridge	6260-H102A
10-56-02004	OPTI-MAX Replacement Cartridges, 1/8" Ball, 2/Pack	n/a
For 880/980 SERIES		
42-56-01544	OPTI-MAX Inlet Housing & Cartridge	6560-H405A
42-56-01546	OPTI-MAX Outlet Housing	

## seals

Optimize P/N	Optimize Description	jasco P/N
For BIP1, PU-2080, 2085, 2089, 1580-c02, 980, 880, 1586, 2086, 986		
42-36-01357	OPTI-SEAL UHMW-PE Piston Seal	6560-H148, H149A
42-36-01358	OPTI-SEAL UHMW-PE Piston Seal, 10/pk	n/a
42-36-01355	ITB PTFE Piston Seal	6560-H148-PE, H149A
42-36-01356	ITB PTFE Piston Seal, 10/Pack	n/a

## pistons

Optimize P/N	Optimize Description	jasco P/N
42-34-01354	Sapphire Piston	2305-0005A

## lamps

Optimize P/N	Optimize Description	JASCO P/N
42-32-02841	DX2 DoubleLife Lamp, D2, 870, 875	5330-0097
42-32-02842	DX2 DoubleLife Lamp, D2, UV-1570, -1575, -2075, 975 B/C Series	5330-0091