


Instrument Replacement Parts

thermo



SURVEYOR L-PUMP

optimax check valves

optimize P/N	Optimize Description	Thermo P/N
50-46-00893	OPTI-MAX Outlet Cartridge Check Valve, 1/16" Ceramic Ball & Seat, SS Cartridge	00950-30021
50-46-01983	OPTI-MAX Inlet Cartridge Check Valve, 1/16" Ruby Ball/Sapphire Seat, SS Cartridge	00950-30026
50-47-00893	OPTI-MAX Outlet Cartridge Check Valve, 1/16" Ruby Ball & Sapphire Seat, SS Cartridge	n/a
50-47-01983	OPTI-MAX Inlet Cartridge Check Valve, 1/16" Ruby Ball & Sapphire Seat, SS Cartridge	n/a

seals

Optimize P/N	Optimize Description	Thermo P/N
50-36-00865	OPTI-SEAL, UHMW-PE Piston Seal	00950-30004
50-36-00866	OPTI-SEAL, UHMW-PE Piston Seal, 10/pk	n/a
50-36-00863	ITB PTFE Piston Seal	00107-18008
50-36-00864	ITB PTFE Piston Seal, 10/pk	n/a

pistons

Optimize P/N	Optimize Description	Thermo P/N
50-34-02650	Sapphire Piston	00950-30020

optimize P/N	Optimize Description	Thermo P/N
50-48-00893	OPTI-MAX Outlet Cartridge Check Valve, 1/16" Ceramic Ball & Seat, PEEK Cartridge	n/a
50-48-01983	OPTI-MAX Inlet Cartridge Check Valve, 1/16" Ceramic Ball & Seat, PEEK Cartridge	n/a
50-49-00893	OPTI-MAX Outlet Cartridge Check Valve, 1/16" Ruby Ball & Sapphire Seat, PEEK Cartridge	n/a
50-49-01983	OPTI-MAX Inlet Cartridge Check Valve, 1/16" Ruby Ball & Sapphire Seat, PEEK Cartridge	n/a



optimize technologies

13993 Fir Street
Oregon City, OR 97045

TF: 800.669.9015 | Tel: 503.557.9994 | Fax: 503.557.9995

www.optimize-tech.com