







DIY EXP[®]2 TI-LOK™ FITTING & ADAPTER KITS

MAKE YOUR OWN HAND-TIGHT UHPLC CONFIGURATIONS

Any pressure. Any tubing. Any length. Any port. Any combination.
The possibilities are endless...

<p>EXP[®]2 TI-LOK™ 10-32 FITTINGS</p>	 <p>For All Tubing</p> <p>For All Tubing</p>	 <p>Any 1/16" OD Tubing</p> <p>EXP[®]2 Driver</p> <p>Any 1/16" OD Tubing</p>
<p>EXP[®]2 TI-LOK™ 6-32 & 6-40 FITTINGS</p>	 <p>For Fused Silica Tubing</p> <p>For All Tubing</p>	 <p>OPTI-PEEK™ Sleeve</p> <p>360µm OD Fused Silica</p> <p>Any 1/32" OD Tubing</p> <p>EXP[®]2 Driver</p>
<p>EXP[®]2 TI-LOK™ 10-32 ADAPTERS</p>	 <p>For Fused Silica Tubing</p> <p>For All Tubing</p>	 <p>360µm OD Fused Silica</p> <p>Any 1/32" OD Tubing</p> <p>EXP[®]2 Driver</p>

For more applications such as sample loops and packed capillary columns, please contact Optimize for more information.